




## CERTIFICATE OF CONFORMITY

We hereby certify that the below listed products are manufactured and inspected according to our approved documents and procedures before they are released for delivery. During inspection the products have to fulfill the set quality criteria and performance specifications.

Product Code:	Product Name:	Lot no:	Exp:
97002620B	KF Flex 24 DW tip comb & plate, purity certified	2220010	2027-07

Our Quality System has been found to conform to the Quality System Standard ISO 9001.

In Joensuu, on 21st July 2022

  
Heidi Salokangas  
Quality Manager  
Laboratory Plastics Essentials  
Thermo Fisher Scientific Oy