

Quality Assurance Certificate

Sartobind Q
Strong basic anion exchanger
MA 100

Order No.: 93IEXQ42BC-12
Lot No.: 240111133

This product is developed, produced and distributed according to a Quality Management System that is certified for compliance with DIN / ISO 9001.
This product has passed Sartorius Stedim Biotech's in-house tests and thus meets Sartorius Stedim Biotech's stringent quality control standards.

This product is tested for binding capacity and flow rate. The membrane is tested for binding capacity, thickness, evenness and flow rate.

Membrane type	:	Quaternary ammonium
Flow rate membrane	:	> 80 ml / (cm ² x min x 0.1 MPa)
Protein binding capacity at 10% breakthrough (Reference protein Bovine serum albumin)	:	> 0.50 mg/cm ²

Note

This product is supplied as non-sterile. The membrane contains glycerol.

2025-01-28
Date



Dr. Anna Vreemann
Site Quality Manager

Manufactured by Sartorius Stedim Biotech GmbH
37070 Goettingen, Germany
Phone: + 49.551.308.0
Fax: + 49.551.308.3289
www.sartorius.com