

## Quality Assurance Certificate

**Minisart® SRP**

**Order no.** 17575-----K  
**Use before** 04 / 2030  
**Pore size** 0.20 µm  
**Lot no.** 250711103

Minisart is manufactured by Sartorius in accordance with established quality standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no. alone or lot no. with one of the following: pore size, ID no., material type, order no. are printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biosafe according to the USP Class VI Biological tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.72 for non-fiber-releasing filters.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

*Burst Pressure Test*

*Bubble Point Test*

*Pressure Hold Test*

*Flow Rate Performance*

*Sterile Filtration Capability (Bacterial Challenge Test)*

2025-04-22

Date



Dr. Matthias Reher  
Site Quality Manager

Manufactured by Sartorius Stedim Biotech GmbH  
37070 Goettingen, Germany

Phone +49.551.308.0  
Fax +49.551.308.3289

[www.sartorius.com](http://www.sartorius.com)