

Sampling bag SteriBag

Item no.	Description
5344-0060	Sampling bag SteriBag, 60ml, w/o lab.area, 127x76mm
5344-0061	Sampling bag SteriBag, 60ml, w/ lab.area, 127x76mm
5344-0180	Sampling bag SteriBag, 150ml, w/o lab.area, 178x76mm
5344-0181	Sampling bag SteriBag, 150ml, w/ lab.area, 178x76mm
5344-0400	Sampling bag SteriBag StandUp, 400 ml, w/ lab.area
5344-0540	Sampling bag SteriBag, 450ml, w/o lab.area, 229x114mm
5344-0541	Sampling bag SteriBag, 450ml, w/ lab.area, 229x114mm
5344-0600	Sampling bag SteriBag StandUp, 600ml, w/o lab.area
5344-0720	Sampling bag SteriBag, 650ml, w/o lab.area, 229x140mm
5344-0721	Sampling bag SteriBag, 650ml, w/ lab.area, 229x140mm
5344-0810	Sampling bag SteriBag, 900ml, w/o lab.area, 305x114mm
5344-0811	Sampling bag SteriBag, 900ml, w/ lab.area, 305x114mm
5344-1650	Sampling bag SteriBag, 1650ml, w/o lab.area, 305x178mm
5344-1651	Sampling bag SteriBag, 1650ml, w/ lab.area, 305x178mm
5344-4000	SteriBag 4000 ml, w/ lab.area
5344-5500	SteriBag 5500 ml, w/o lab.area
5344-7008	Sampling bag SteriBag blue, 650ml, 229x140mm
5344-7012	Sampling bag SteriBag blue, 1650ml, 305x178mm
5344-8010	Sampling bag SteriBag cleanroom, 800ml, 305x114mm
5344-8012	Sampling bag SteriBag cleanroom, 1650ml, 305x178mm
9990-0354	SteriBag 25,4x45,7cm, l.a sterile

Our supplier certifies random tests are carried out for aerobic and anaerobic bacteria as well as yeast and molds. Sterility tests, based on USP 39 Microbiological Tests / <71> Sterility Tests, are performed to guaranty sterility of the process.

These tests confirm process sterility during the production of following LOT no.:

LOT:	Expiry date:	LOT:	Expiry date:
J104646	2030-05	J104654	2030-05
J104647	2030-05	J104655	2030-05
J104649	2030-05	J104656	2030-05
J104650	2030-05	J104657	2030-05
J104651	2030-05	J104658	2030-05
J104652	2030-05	J104661	2030-05
J104653	2030-05	J104826	2030-05

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

16. June 2025

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director

