

Corning Incorporated - Life Sciences
 271 County Route 64
 Big Flats NY 14814 USA
www.corning.com/lifesciences
 Refer to website for regional contact information.

Page: 1 / 1

| | | | |
|------------------------|---|-------------------------|--------------|
| Product Name | : CoolRack XT 5ml holds 12 x 5ml centrifuge tubes | | |
| Catalog Number | : 432060 | Manufacture Date | : 2016-05-24 |
| Lot ID | : 04216027 | | |
| Expiration Date | : Unconfirmed Validity | | |

Quality Management System - Complies with the current version of the ISO 9001 Standard.

Quality Control Testing - Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release.

- This product met Corning Incorporated - Life Sciences' high standards of quality at the time of batch/lot release.



David R Patrick
 Quality Manager