



Fisher Scientific Company LLC.  
300 Industry Drive Pittsburgh  
Pittsburgh, PA 15275-1014  
Tel: 724-517-2400  
Fax: 724-517-1546  
fishersci.com

Catalog Number: 0666247

Lot Number: SAMPLE

Description: Fisherbrand™ Traceable™ Two-Channel Benchtop Timer with Dual-Line LCD

### **Certificate of Conformance**

Our supplier hereby certifies that the items advertised in our catalogs and/or website conform to the specifications listed. Any changes made to the specifications after our catalog print date are on file and available upon request or you may visit [www.coleparmer.com](http://www.coleparmer.com) for the most current data.

Our supplier cannot and does not assume any responsibility for any failure of the product to meet this Certificate of Conformance that results from any handling or use subsequent to opening the package. Except as otherwise set forth in writing by our Supplier, there are no other warranties, expressed or implied, that apply to this certificate or the subject product(s).

Our supplier does not test or analyze these products for any specified regulatory requirements. Therefore, customers and end-users must make their own determinations ensuring the uses of these products are safe, lawful and suitable for their intended applications.

This Certificate has been compiled with data that has been provided by the manufacturer of the product. The Fisher Scientific channel has not independently verified such data. This document has been produced electronically and is valid without a signature.