

Certificate of Irradiation

Product	Nunc TripleFlasks Vent/Close, 500 cm ²	Date 2025-09-02
		Page 1
Item Number	132867	
Lot Number	186585	
Specified Dose Range	22.2 - 40.3 kGy	
SAL	10E-6	
Quantity	200 boxes	

Thermo Fisher Scientific confirms that the product supplied has been inspected and tested in accordance to relevant specifications.

Irradiation lot no.	PalletID	Date of irradiation	Dose kGy
EA__1724219A	DK00506171	2025-07-13	28.5
EA__1723960B	DK00506172	2025-07-13	29.7
EA__1724003B	DK00506173	2025-07-13	27.9
EA__1723987B	DK00506174	2025-07-13	27.9
EA__1723505B	DK00506175	2025-07-23	27.8
EA__1724345A	DK00506176	2025-08-08	27.3
EA__1724367B	DK00506177	2025-08-08	27.4
EA__1724322A	DK00506178	2025-08-08	28.5
EA__1724344A	DK00506179	2025-08-08	27.3
EC__1300378A	DK00506180	2025-08-25	26.6
EC__1300378A	DK00506180	2025-08-25	26.6

Sterility is obtained through irradiation according to ISO 11137 "Sterilization of health care products - Radiation".

Quality assurance of the product is performed in accordance with the requirements of the Quality Management System ISO 9001:2015 "Quality Management System - Requirements" and ISO 13485:2016 "Medical Devices - Quality Management Systems - System requirements for regulatory purposes".

Roskilde Site



Per Bernack
Senior Quality Assurance Manager

Please note that this is an electronically generated signature