

Certificate of Irradiation

Product Name: Reagent Reservoir
Description: Reagent Reservoir, 10mL, Polystyrene, Sterile, 1 per Bag, Pack of 100
Heathrow Part Number: 120639
Lot Number: 601543
Expiration Date: 01/31/2029

Sterilization

Irradiation Lot Number: 00008123
Irradiation Finished Date: 12/30/2024
Dose Inspection Result: 25.02 kGy
Specified Dose: 25-50 kGy

Product has been irradiated with high-energy electron beam.

Nicole Smokowicz

Quality Assurance Officer