

## Certificate of Irradiation

**Product Name:** Reagent Reservoir  
**Description:** Reagent Reservoir, 25mL, Polystyrene, Sterile, 1 per Bag, Pack of 100  
**Heathrow Part Number:** 120641  
**Lot Number:** 601543  
**Expiration Date:** 01/31/2029

### **Sterilization**

Irradiation Lot Number: 00008040  
Irradiation Finished Date: 12/23/2024  
Dose Inspection Result: 25.02 kGy  
Specified Dose: 25-50 kGy

Product has been irradiated with high-energy electron beam.

*Nicole Smokowicz*

Quality Assurance Officer