

Certificate of Conformance and Analysis

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

Product Code:	PNHS-66
Lot Number:	022330
Manufacturing Date:	2024.09.23
Quantity of Particles Po:	9.11 millions/m ²
Quantity of Fibers Po:	0.07 thousands/m ²
Sorbency Rate:	1.0 seconds
Non-volatile Residue in IPA:	0.015 g/m ²
Non-volatile Residue in DI:	0.006 g/m ²

Quality Department
Contec Cleanroom Technology (Suzhou) Co., Ltd.
2024-10-16



NORTH AMERICA

Contec, Inc.
525 Locust Grove
Spartanburg, SC 29303 USA
tel: +1 864 503 8333
fax: +1 864 503 8444

EUROPE

Contec Europe
ZI Du Prat, Avenue Paul Duplex
56000 VANNES France
tel: +33 (0)2 97 43 76 98
fax: +33 (0) 97 43 76 86

ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.
17 Longyun Road
Suzhou Industrial Park, Suzhou 215024 China
tel: +86 512 6274 4050
fax: +86 512 6274 4051

web: www.contecinc.com
e-mail: info@contecinc.com