



## PRODUCT CERTIFICATE

Product: Contec NeutraKlean  
Product Code: FBC103NK  
Product Description: Filtered Neutral Detergent in purified water (EP) 1L Capped  
Batch Number: 250503083  
Manufacture Date: May / 2025  
Expiry Date: May / 2027

### ANALYSIS

Test	Specification	Results
Colour:	Colourless	Complies
Clarity:	Clear	Complies
Filtration:	Filtered to 0.2 microns	Complies
SG at 20°C:	0.995 - 1.010	1.000
pH:	6.0 to 8.0	6.68

Manufactured product via Quality Systems certified to ISO 9001:2015, tested in accordance with documented quality procedures and approved when required specifications are met.

Name: 1: Nicola Hunter 2: Kobey Robson  
Position: 1: Quality Assurance Officer 2: Quality Assurance Officer  
Date: 1: 18 JUL 25 2: 21 JUL 25  
Authorised Signature: 1:  2: 

For and on behalf of Contec Inc

COA55 Rev 3

**Manufactured by:**  
Contec Cleanroom (UK) Ltd  
Unit 6A Wansbeck Business Park  
Ashington  
UK

**America**  
Contec Inc  
P.O.Box 530  
Spartanburg SC  
USA

**Europe**  
Contec Inc  
ZI du Prat RP 3707  
56037 VANNES  
France

**China**  
Contec Cleanroom Technology (Suzhou) Co. Ltd  
No. 17 Longyun Road  
Suzhou 215024  
China

[www.contecinc.com](http://www.contecinc.com)  
[infoeu@contecinc.com](mailto:infoeu@contecinc.com)