



## Certificate of Conformance and Analysis

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

---

Product Code:	492222-813
Lot Number:	1079949-0-0-1

---

### Analysis

Biaxial Shake Fibers:	0.247 thousands/m <sup>2</sup>
Biaxial Shake Particles 0.5µM:	10.9 millions/m <sup>2</sup>
Sorbency Rate:	1.0 seconds
Nonvolatile residue in DI:	0.0155 g/m <sup>2</sup>
Nonvolatile residue in IPA:	0.031 g/m <sup>2</sup>

---

A handwritten signature in black ink, appearing to read "C. Duffee".

QC Supervisor  
coc@contecinc.com

Tuesday, July 23, 2024

#### NORTH AMERICA

Contec, Inc.  
525 Locust Grove  
Spartanburg, SC 29303 USA  
tel: +1 864 503 8333  
fax: +1 864 503 8444

#### EUROPE

Contec Europe  
ZI Du Prat, Avenue Paul Duplex  
56000 VANNES France  
tel: +33 (0)2 97 43 76 98  
fax: +33 (0) 97 43 76 86

#### ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.  
17 Longyun Road  
Suzhou Industrial Park, Suzhou 215024 China  
tel: +86 512 6274 4050  
fax: +86 512 6274 4051

Page 1 of 1

web: [www.contecinc.com](http://www.contecinc.com)  
e-mail: [info@contecinc.com](mailto:info@contecinc.com)

DID: 9361B