

PRODUCT CERTIFICATE

Product: Contec ProChlor

Product Code: FBC102PC

Product Description: Filtered Stabilised Hypochlorous Acid in purified water (EP) 1L Capped

Batch Number: 240602796

Manufacture Date: JUN / 2024

Expiry Date: JUN / 2026

ANALYSIS



Test	Specification	Results
Colour:	Colourless	Complies
Clarity:	Clear	Complies
Filtration:	Filtered to 0.2 microns	Complies
SG at 20°C:	0.990 – 1.010	1.000
Available chlorine:	>1000ppm	2184
pH at 20°C :	3.0 – 6.0	3.52

Manufactured product via Quality Systems certified to ISO 9001:2015, tested in accordance with documented quality procedures and approved when required specifications are met.

Name: 1: Kobey Robson 2: Nicola Hunter

Position: 1: Quality Assurance Officer 2: Quality Assurance Officer

Date: 1: 25 JUL 24 2: 25 JUL 24

Authorised Signature: 1:  2: 

For and on behalf of Contec Inc

COA51 Rev 2

Manufactured by:
Contec Cleanroom (UK) Ltd
Unit 6A Wansbeck Business Park
Ashington
UK

America
Contec Inc
P.O.Box 530
Spartanburg SC
USA

Europe
Contec Inc
ZI du Prat RP 3707
56037 VANNES
France

China
Contec Cleanroom Technology (Suzhou) Co. Ltd
No. 17 Longyun Road
Suzhou 215024
China

www.contecinc.com
infoeu@contecinc.com