



Certificate of Conformance and Analysis

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

Product Code:	495353-903
Lot Number:	1077781-0-0-1
Manufacturing Date:	07/10/2024

Analysis

Biaxial Shake Fibers:	0.095 thousands/m ²
Biaxial Shake Particles 0.5µM:	7.6 millions/m ²
Sorbency Rate:	1.0 seconds
Nonvolatile residue in DI:	0.00 g/m ²
Nonvolatile residue in IPA:	0.0219 g/m ²

A handwritten signature in black ink, appearing to read "C. Duffee".

QC Supervisor
coc@contecinc.com

Saturday, July 13, 2024

NORTH AMERICA

Contec, Inc.
525 Locust Grove
Spartanburg, SC 29303 USA
tel: +1 864 503 8333
fax: +1 864 503 8444

EUROPE

Contec Europe
ZI Du Prat, Avenue Paul Dupleix
56000 VANNES France
tel: +33 (0)2 97 43 76 98
fax: +33 (0) 97 43 76 86

ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.
17 Longyun Road
Suzhou Industrial Park, Suzhou 215024 China
tel: +86 512 6274 4050
fax: +86 512 6274 4051

web: www.contecinc.com
e-mail: info@contecinc.com

Page 1 of 1

DID: 8121B