

Certificate of Conformity

Konformitätsbescheinigung
Certificat de Conformité
Certificato di conformità
Certificado de conformidad
Certificado de Conformidade



Page 1 of 1

Part Number	M1040B - 100
Description	Metal Detectable Long Handled Scoop - Blue
Batch Number	065817
Date of Manufacture	9 April 2019
Expiry Date	9 April 2024
Material	PP Polypropylene (with Metal Detectable Filler)
Certificate Issued	02 May 2019

- We confirm that the goods are manufactured with Polypropylene that conforms to:
 - PP conforms to FDA 21 CFR 177.1520
 - Blue colourant conforms to FDA 21 CFR 178.3297
 - Metal Detectable Filler conforms to FDA 21 CFR 186.1374
 - Conforms to EU Regulations 10/2011
 - Conforms to EC Regulations 1935/2004
- An FDA acceptable blue dye (less than 1%) is added during manufacture to make the product blue.
- The ingredients used in the manufacture of this product are certified BSE/TSE free. During manufacture this item does not come into contact with any products or processes that may pose a BSE/TSE risk.
- This product is moulded in a class 100,000 Cleanroom.
- This product is packed in a class 10,000 Cleanroom
- This product is individually bagged.
- The information detailed above has been given to us by our suppliers, it is passed on in good faith and is to the best of our knowledge correct.

Signed on behalf of Sampling Systems Ltd.

A handwritten signature in blue ink, appearing to read 'Kuldeep Flora', written over a light blue horizontal line.

Kuldeep Flora

