

Certificate of Quality

Item Description: OvermoldMillex13 0.20mPhobicPTFEs100pkR

Catalogue Number: SLFGX13NL

Lot Number: 0000360749

ISO 9001 Quality Standard

EMD Millipore Corporation certifies that the above-referenced product:

Was manufactured (or supplied) by a facility whose Quality Management System is approved by an accredited registering body to the appropriate ISO 9001 Quality Systems Standard.

Was manufactured in accordance with applicable document standard procedures.

Catalogue Number	Pack Size	Membrane Type	Pore Size (µm)	Diameter	Typical Methanol Flow Rate at 21 degrees C and 2.1 bar (30 psi)	Bubble point
SLFGX13NL	100	Millex®-FG hydrophobic PTFE membrane	0.2	13mm	53 mL/min	≥ 0.8 bar (≥ 11 psi)
SLFGX13TL	100					
SLFGX13NK	1000					
SLFG033NS	50					
SLFG033NB	250	Millex®-FH hydrophobic PTFE membrane	0.2	33mm	259 mL/min	≥ 0.8 bar (≥ 12 psi)
SLFG033NK	1000					
SLFHX13NL	100		0.45	13mm	137 mL/min	≥ 0.6 bar (≥ 8 psi)
SLFHX13TL	100					
SLFHX13NK	1000					
SLFH033NS	50					
SLFH033NB	250			33mm	407 mL/min	≥ 0.6 bar (≥ 8 psi)
SLFH033NK	1000					

The vibrant M, Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

