



CERTIFICATE OF CONFORMANCE

Product Part Number: 825-12AM
Lot Number: 050571-000
Irradiation Reference Number: N/A
Lot Quantity: 96 EACH
Date of manufacture: 07/10/2025
Shelf Life: N/A
Description: CLEANROOM POLYIMIDE TAPE (24/CS)
Size -
Width: ½"
Length: 36 YARDS
Height: N/A

This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.

8/14/2025

Authorized Signature
Quality Control Department

Date

CC: With Product
Customer File
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010 Tel. (310) 784-6990 • FAX (310) 784-6980