



CERTIFICATE OF CONFORMANCE

Product Part Number: CR100PC-1WH
Lot Number: 049665-000
Irradiation Reference Number: N/A
Lot Quantity: 720 EACH
Date of manufacture: 06/17/2025
Shelf Life: N/A
Description: CLEANROOM TAPE, WHITE, VINYL (36/CS)
Size -
Width: 1"
Length: 36 YARDS
Height: N/A

This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.

7/17/2025

Authorized Signature
Quality Control Department

Date

CC: With Product
Customer File
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010 Tel. (310) 784-6990 • FAX (310) 784-6980