



CERTIFICATE OF CONFORMANCE

Product Part Number: 825-1AM
Lot Number: 049116-000
Irradiation Reference Number: N/A
Lot Quantity: 540 EACH
Date of manufacture: 04/29/2025
Shelf Life: N/A
Description: CLEANROOM POLYIMIDE TAPE (36/CS)
Size -
Width: 1"
Length: 36 YARDS
Height: N/A

This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.

5/15/2025

Authorized Signature
Quality Control Department

Date

CC: With Product
Customer File
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010 Tel. (310) 784-6990 • FAX (310) 784-6980