

CERTIFICATE OF CONFORMANCE

Product Part Number: 825-1AM1
Lot Number: 023805-000
Irradiation Reference Number: N/A
Lot Quantity: 36 EACH
Date of manufacture: 02/02/2021
Shelf Life: N/A
Description: CLEANROOM POLYIMIDE TAPE (36/CS)
Size -
Width: 1"
Length: 22 YARDS
Height: N/A

This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.

Sabito

02/25/2021

Authorized Signature
Quality Control Department

Date

CC: With Product
Customer File
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010 Tel. (310) 784-6990 • FAX (310) 784-6980