

CERTIFICATE OF ANALYSIS

Thermo Fisher Scientific's Quality System has been found to conform to Quality Management System Standard ISO9001:2015 by Intertek Global Certificate Number. CERT-0120633

Catalogue Number	BP665
Lot Number	254567
Description	PHOSPHATE BUFFER SALINE, 10X READY CONCENTRATE
CAS Number	7647-14-5
Quality Test/Release Date	23/Oct/2025
Suggested retest date	31/Mar/2028
Country of Origin	United States
Declaration of Origin	Inorganic
BSE/TSE	No animal products are used as starting raw material ingredients, or used in processing, including lubricants, processing aids, or any other material that might migrate to the finished product.

Result Name	Units	Specifications	Test Value
APPEARANCE		REPORT	White granular powder
CHLORIDE	MOLES/LITER	Inclusive Between 0.136 - 0.144	0.140
CONDUCTIVITY AT 25 DEG C	MICROMHOS/CM	Inclusive Between 14000 - 17800	16329
DNASE	PASS/FAIL	= NOT DETECTED	NOT DETECTED
PROTEASE	PASS/FAIL	= NOT DETECTED	NOT DETECTED
RNASE	PASS/FAIL	= NOT DETECTED	NOT DETECTED
SOLUBILITY	PASS/FAIL	= PASS TEST	PASS TEST
WEIGHT OF SOLID IN PACKAGE	GRAMS	Inclusive Between 97.3 - 100.3	100.1
pH OF 10X SOLUTION		Inclusive Between 6.7 - 6.9	6.9
pH OF A 1X SOLUTION		Inclusive Between 7.3 - 7.5	7.4



Ash Frey
Supervisor, Quality Control

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third-party data or information associated with the product. Products are for research use or further manufacturing. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

Note: The data listed is valid for all package sizes of this lot of this product, expressed as an extension of the catalogue number listed above.